

GE Healthcare

Aestiva/5 anesthesia machine

More than superior ventilation

Features

Superior ventilation: 7900 SmartVent™

- Volume Mode, pressure
- Volume Mode, Pressure Control Mode, Pressure Support (PSVPro®), Synchronized Intermittent Mandatory Ventilation (SIMV), electronic PEEP
- Tidal volume compensation
- One motion from mechanical to manual mode
- Two key presses to total standby; end case
- Cardiac bypass case mode

Open systems architecture

- Lower overall height
- User configurable drawers/shelving

Innovative patient breathing system

- Eight machine hoses/cables integrated
- "No tools" disassembly of components
- Autoclavable and latex-free
- Responsive location of common gas outlet



Aestiva®/5
Two vaporizer configuration



Aestiva/5
Three vaporizer configuration

Improved low flow/reduced life cycle costs

- Fresh gas flow compensation—automatically
- Smooth, faster acting fresh gas flow control
- Minimum O₂ flow of 50 mL
- Dual air flow tube for low flow
- Two scheduled maintenance checks per year



Physical Specifications

Dimensions

	2 vaporizer configuration	3 vaporizer configuration
Height:	135.8 cm/53.4 in	135.8 cm/53.4 in
Width:	75 cm/29.5 in	93 cm/36.6 in
Depth:	83 cm/32.7 in	83 cm/32.7 in
Weight:	Approximately 136 kg/300 lb	Approximately 154 kg/340 lb

Top shelves (optional)

	2 vaporizer configuration	3 vaporizer configuration
Weight limit:	46 kg/100 lb	46 kg/100 lb
Width:	47.5, 67.5 or 87.5 cm/ 18.7, 26.6 or 34.4 in	87.5 or 67.5 cm/ 34.4 or 26.6 in
Depth:	41 cm/16.1 in	41 cm/16.1 in

Work surface

Height:	87.6 cm/34.5 in
Width:	47 cm/18.5 in
Depth:	31.5 cm/12.4 in

Folding side shelf (optional)

Height:	87.5 cm/34.5 in
Width:	26.5 cm/10.4 in
Depth:	31.5 cm/12.4 in
Weight limit:	11.3 kg/25 lb

DIN rail (optional)

Side of tabletop:	30 cm/12 in
Side of machine:	23.5 cm/9.25 in

Top drawer (1 standard)—locking (internal dimensions)


Height:	10.5 cm/4.1 in
Width:	38.5 cm/15.2 in
Depth:	26 cm/10.2 in

Lower drawers (optional)*

Height:	14.5 cm/5.7 in
Width:	38.5 cm/15.2 in
Depth:	26 cm/10.2 in

Lower shelves (optional)*

Heights:	9.2 cm/3.7 in	13.2 cm/5.2 in
	20.6 cm/8.2 in	24.6 cm/9.8 in
	28.6 cm/11.4 in	36 cm/14.4 in
Width:	42.5 cm/16.75 in	42.5 cm/16.75 in
Depth:	36 cm/14 in	36 cm/14 in




Space for additional shelves and drawers

* Lower cabinet can be configured with a variety of shelf and drawer combinations

Absorber arms		
	Adjustable	Non-adjustable
Arm length:	30.5 cm/12 in	25.4 cm/10 in
Bag arm height:	87 cm/34.3 in 104 cm/40.9 in	91.5 cm/36 in
Absorber rotation:	85°	85°
Ventilator screen		
Height:	7.6 cm/3 in	
Width:	15.2 cm/6 in	
Casters		
Diameter:	12.5 cm/5 in	
Brakes:	Single foot lever locks and unlocks two front casters	

Ventilator operating specifications



Ventilation operating modes	
Volume Control	
Pressure Control	
Synchronized Intermittent Mandatory Ventilation (SIMV)	
Pressure Support (PSVPro) with Apnea Backup ventilation — (optional)	

Ventilator (V_T) parameter ranges	
Tidal volume range:	20 to 1500 mL (Volume Control and SIMV modes) 5 to 1500 mL (Pressure Control Mode)
Incremental settings:	20 to 100 mL (increments of 5 mL) 100 to 300 mL (increments of 10 mL) 300 to 1000 mL (increments of 25 mL) 1000 to 1500 mL (increments of 50 mL)
Minute volume range:	0 to 99.9 L/min
Pressure ($P_{inspired}$) range:	5 to 60 cm H ₂ O (increments of 1 cm H ₂ O)
Pressure (P_{lim}) range:	12 to 100 cm H ₂ O (increments of 1 cm H ₂ O)
Pressure ($P_{support}$) range:	Off, 2 to 40 cm H ₂ O (increments of 1 cm H ₂ O)
Rate:	4 to 100 breaths per minute for Volume Control and Pressure Control; 2 to 60 breaths per minute for SIMV, PSVPro and SIMV-PC+PSV (increments of 1 breath per minute)
Inspiratory/expiratory ratio:	2:1 to 1:8 (increments of 0.5)
Inspiratory time:	0.2 to 5.0 seconds (increments of 0.1 seconds) (SIMV and PSV Pro)
Trigger window:	0 to 80% (increments of 5%)
Flow trigger:	0.2 to 1.0 L/min (increments of 0.2 L/min) 1 to 10 L/min (increments of 0.5 L/min)
Inspiration termination level:	5 to 75% (increments of 5%)
Backup mode delay:	10 to 30 seconds (increments of 5 seconds)

Positive End Expiratory Pressure (PEEP)		Alarm settings	
Type:	Integrated, electronically controlled	Tidal volume (V_T):	Low: OFF, 0 to 1500 mL. High: 20 to 1600 mL, OFF
Range:	OFF, 4 to 30 cm H ₂ O (increments of 1 cm H ₂ O)	Minute volume (V_E):	Low: OFF, 0 to 10 L/min High: 0 to 30 L/min, OFF
Ventilator performance		Inspired oxygen (FiO ₂):	Low: 18 to 100% High: 18 to 100%, OFF
Pressure range at inlet:	240 kPa to 700 kPa/ 35 psig to 100 psig	Apnea alarm:	<i>Mechanical ventilation ON:</i> < 5 mL breath measured in 30 seconds <i>Mechanical ventilation OFF:</i> < 5 mL breath measured in 30 seconds
Peak gas flow:	120 L/min + fresh gas flow	Low airway pressure:	4 cm H ₂ O above PEEP
Flow valve range:	1 to 120 L/min	High pressure:	12 to 100 cm H ₂ O (increments of 1 cm H ₂ O)
Flow compensation range:	200 mL/min to 15 L/min	Sustained airway pressure:	<i>Mechanical ventilation ON:</i> $P_{\text{limit}} < 30$ cm H ₂ O, the sustained limit is 6 cm H ₂ O $P_{\text{limit}} 30$ to 60 cm H ₂ O, the sustained limit is 20% of P_{limit} $P_{\text{limit}} > 60$ cm H ₂ O, the sustained limit is 12 cm H ₂ O <i>PEEP and mechanical ventilation ON:</i> Sustained limit increases by PEEP minus 2 cm H ₂ O <i>Mechanical ventilation OFF:</i> $P_{\text{limit}} < 60$ cm H ₂ O, the sustained limit is 50% of P_{limit} $P_{\text{limit}} > 60$ cm H ₂ O, the sustained limit is 30 cm H ₂ O
Ventilator monitoring		Subatmospheric pressure:	Paw < -10 cm H ₂ O
Expiratory minute volume range:	0 to 99.9 L/min	Alarm silence countdown timer:	120 to 0 seconds
Expiratory tidal volume range:	0 to 1500 mL		
O ₂ %:	5 to 110%		
Peak pressure:	-20 to 120 cm H ₂ O		
Mean pressure:	-20 to 120 cm H ₂ O		
Plateau pressure:	0 to 120 cm H ₂ O		
Pressure waveform sweep speed:	4 to 25 breaths per minute (0 to 15 seconds) 26 to 75 breaths per minute (0 to 5 seconds) 75 breaths per minute (0 to 3 seconds)		
Ventilator accuracy			
Delivery/monitoring accuracy			
Volume delivery:	> 210 mL = better than 7% < 210 mL = better than 15 mL < 60 mL = better than 10 mL.		
Pressure delivery:	±10% or ±3 cm H ₂ O		
PEEP delivery:	±1.5 cm H ₂ O		
Volume monitoring:	> 210 mL = better than 9% < 210 mL = better than 18 mL < 60 mL = better than 10 mL.		
Pressure monitoring:	±5% or ±2 cm H ₂ O		

Ventilator components

Flow transducer

Type:	Variable orifice flow sensor
Dimensions:	22 mm OD and 15 mm ID
Location:	Inspiratory outlet and expiratory inlet


Optional autoclavable sensor available

Oxygen sensor


Type:	Galvanic fuel cell
Life cycle:	Approximately 18 months (dependent on usage)

Anesthetic agent delivery

Vaporizers:	Tec 4, Tec 5, Tec 6 Plus, Tec 7
Number of positions:	2 or 3
Mounting:	Tool-free installation Selectatec® manifold interlocks and isolates vaporizers



Tec 7



Tec 6 Plus

Electrical specifications

Current leakage

120 V:	< 300µA
--------	---------

Light package

Task light:	12 V, 3 lamps, type 194, .270A each
Goose neck (optional):	12 V, type 1815, .200A

Power and battery backup

Power input:	120 Vac, 60 Hz, 10A
Backup power:	Demonstrated battery backup time under typical operating conditions is 45 minutes when fully charged
Battery type:	Internal rechargeable sealed lead acid
Power cord:	Length: 5 m/16.4 ft Rating: 15A @ 120 Vac

Communication port

Serial interface:	Isolated RS-232C compatible port
-------------------	----------------------------------

Inlet/outlet modules (120 V)

System circuit breakers:	No outlets 5A w/outlets 10A
Outlets (optional):	4 outlets on back, 3-2A, 1-3A individual breakers and 1-5A combined outlet breaker; optional isolation transformer
Auxiliary outlet box (optional):	5 NEMA outlets on dovetail-mounted box, 5-2A breakers, isolation transformer
Tec 6 Plus outlet:	1 IEC 320 located above vaporizer backbar

Pneumatic specifications		Flowmeters	
Internal common gas outlet		O ₂ ranges:	Two tubes: 0.05 to 0.95 L/min and 1 to 15 L/min
Connector:	ISO 22 mm OD and 15 mm ID		Minimum O ₂ flow: 50 mL/min ±25 mL
Auxiliary common gas outlet (optional)		N ₂ O ranges:	Two tubes: 0 to 0.95 L/min and 1 to 10 L/min
Connector:	ISO 22 mm OD and 15 mm ID	Air range:	One tube option: 1 to 15 L/min Two tube option: 0 to 0.95 and 1 to 15 L/min (low flow tube optional)
Gas supply		CO ₂ (optional):	One tube: 0 to 0.5 L/min
Pipeline input range:	240 kPa to 600 kPa/ 35 psig to 88 psig	Heliox range (optional):	One tube: 0 to 15 L/min
Pipeline connections:	DISS-male All fittings available for O ₂ , N ₂ O, and Air, and contain pipeline filter and check valve.	Calibration:	Percent of full scale flow Accuracy (% of flowrate)
Cylinder input:	Pin indexed in accordance with CGA-V-1; contains input filter and check valve Note: Maximum 5 cylinders total; one oxygen required.		100 ±2.5%
Primary regulator diaphragm minimum burst pressure:	2758 kPa/400 psig		90 ±2.5%
Primary regulator nominal output:	< 338 kPa/49 psig Pin indexed cylinder connections		80 ±2.6%
			70 ±2.7%
			60 ±2.9%
			50 ±3.1%
			40 ±3.4%
			30 ±4.0%
			20 ±5.0%
			10 ±8.1%
Gas power outlet (optional)		Calibration conditions:*	20°C/68°F 101.3 kPa/760 mmHg
Connector:	DISS indexed in accordance with CGA-V-5		
Gas:	Oxygen		
Pressure and flow characteristics:	Varies with source		
O₂ controls		Hypoxic guard system	
Method:	Proportionate decrease of N ₂ O, CO ₂ , O ₂ /He with reduction in O ₂ pressure	Type:	Mechanical Link-25™
Supply failure alarm:	Range: 193 kPa to 221 kPa/ 28 psig to 32 psig Sounds at maximum volume every 10 seconds	Range:	Provides a nominal 25% concentration of oxygen in any O ₂ /N ₂ O mixture
O ₂ flush:	Range: 35 to 50 L/min	Materials	
		All materials in contact with patient gas are free of natural rubber latex.	

Environmental specifications

System operation

Temperature:	10° to 40°C/50° to 104°F
Humidity:	15 to 95% relative humidity (non-condensing)
Altitude:	-440 to 3565 m/500 to 800 mmHg

System storage

Temperature:	-25° to 65°C/-13° to 149°F
Humidity:	10 to 100% relative humidity (including condensing)
Altitude:	-440 to 5860 m/375 to 800 mmHg
Oxygen cell storage:	-15° to 50°C/5° to 122°F 10 to 95% relative humidity 500 to 800 mmHg

Electromagnetic compatibility

Immunity:	Complies with all requirements of EN 60601-1-2
Emissions:	CISPR 11 group 1 class B
Approvals:	UL 2601-1, CSA C22.2 #601.1 IEC 601-1 EN 60601-1