

Large Volume Containers with Sliding, Sealing Gasket, and Hinged Lids



Features:

- Adjustable Rotor Opening and Hinged Lids accommodate a variety of sizes of sharps and provide temporary and permanent closure.
- Optional gasketed lids provide leak resistant containment of contents. Containers with gaskets have DOT approved UN3291 for transport of regulated medical waste except for cultures and stocks.
- Containers lock for final disposal.
- Nestable containers save valuable storage space.
- Lids can provide temporary closure when container is not in use to reduce exposure to contents.
- Sliding gasket lids have a round 6" opening to accommodate up to 1000ml bottles.

Code	Description	Dimensions	Ship Case
8980AC	8 Gal Red Container with Clear, Hinged Lid	17.75"H x 11"D x 15.5"W	3/Case
8997	8 Gallon, Red BIOMAX Container with Hinged, Sealing Gasket Lid	17.5"H x 11"D x 15.5"W	10/Case
8980	8 gallon Red Container with Hinged lid	17.5"H x 11"D x 15.5"W	10/Case
8932	12 Gallon, Red BIOMAX Container with Hinged, Sealing Gasket Lid	18.75"H x 12.75"D x 18.25"W	10/Case
8933	12 gallon Red Container with Hinged lid	18.75"H x 13"D x 18"W	10/Case
8991	18 gallon Red Container with Hinged lid	26"H x 12.75"D x 18.25"W	5/Case
8998	18 Gallon, Red BIOMAX Container with Hinged, Sealing Gasket Lid	26"H x 12.75"D x 18.25"W	5/Case
8930SA	30 Gallon, Red Perfusion Waste Container with Hinged, Sealing Gasket Lid	28"H x 22"W x 4.5"D	3/Case
8992H	Floor bracket with wheels for 8, 12 & 18 gallon Large Volume and CHEMOSAFETY containers	22"H x 12.75"D x 22.5"W	1