

# Shuttle Bariatric | HEAVY DUTY

- Elevating legrest models feature a notched, stainless steel ratchet bar to lock the legrests securely in place
- Removable padded armrests and legrests provide extra comfort and make exiting easy
- Double cross-brace for maximum support and durability
- Hemi-height adjustable with a triple axle design
- Attractive gray hammertone finish
- Comfortable nylon upholstery



MDS809750

## ACCESSORIES

### MDS851965SH

Anti-Fold Device

### MDS85197SH

Overhead Anti-Theft Device

### MDS85189SH

Anti-Tip Device

### WCA806965SH

Swing Away Footrest

### WCA806991

Tinnerman Legrest Locks

CHAIR	ARM STYLE	LEG STYLE	SEAT W X D INCHES (CM)	HCPCS	WT. CAP. LB (KG)
<del>MDS809850</del>	Desk Length Removable	Detachable Elevating Legrest	24 x 20 (61 x 51)	K0007+E2202 K0195 or E0990	700 (318)
<b>MDS809750</b>	Desk Length Removable	Detachable Elevating Legrest	26 x 20 (66 x 51)	K0007+E2202 K0195 or E0990	700 (318)
MDS809850	Desk Length Removable	Detachable Elevating Legrest	28 x 20 (71 x 51)	K0007+K0108 K0195 or E0990	700 (318)
MDS809585	Full Length Removable	Detachable Elevating Legrest	30 x 20 (76 x 51)	—	700 (318)

11427