



S-Scort® new-Sentinel

- Small - only 17" x 9" x 5 1/4"
- Light - just 8.5 lbs.
- Powerful - 50mmHg to > 550mmHg (+/- 10%)

This small, light suction unit delivers 30+ lpm clinical airflow at 550+ mmHg. The sealed lead-acid battery and hospital-grade, dual-mode charging system keep the unit constantly ready for use. Variable regulator for infant, child, and adult. Running time is approximately 30-45 minutes.

Dimensions: 17" x 9" x 5 1/4"

Weight: 9 lbs.

Electrical: 120V/60Hz

Product No.: AE-6905

Accessories

Disposable Suction Canister, 1200cc.

Weight: 1 lbs.
Product No.: AE-6974

Dual-Mode, Hospital Grade Charger

Weight: 2 lbs.
Product No.: 80525

Patient Connecting Tubing

Weight: 1 lbs.
Product No.: 43277

Large Bore Suction Tip
Order 5 boxes for \$82.50.

Weight: 2 lbs.
Product No.: AE-6995