

# Stryker Prime™ Series



The Big Wheel can reduce start-up force up to 50 percent and steering effort up to 60 percent.

## Mobile. Safe. Efficient.

Extraordinarily mobile and user-friendly, the Stryker Prime Series stretchers are designed to greatly reduce the physical strain on clinicians while enhancing patient care and comfort. Advanced mobility options such as the Zoom Motorized Drive System or Big Wheel make transporting heavy patients easy, while electric functionality allows patients to adjust their own position without calling a caregiver for assistance. Fully equipped and highly configurable, the Stryker Prime Series allows you to create the ideal solution to meet your specific needs.



Zoom Motorized Drive System on the Stryker Prime stretcher virtually eliminates manual pushing.



Lift Assist™ Backrest uses a drop-seat design which uses the patient's weight to help raise the head of the bed while cradling the patient to help reduce boosting.

## Medical

**Advanced mobility, with optional Big Wheel® or Zoom® Motorized Drive System, helps make transport safe and efficient**

**NEW electric-powered base and litter option available on all mobility solutions**

**Industry's first stretcher exit system, alerts the nurse to a potential patient fall**

**Optional integrated scale system provides accurate, repeatable weighings**

**Pioneer® mattress offers superior pressure redistribution and comfort**

**Power-washable design with 700 lb capacity**



Electric functionality, including exclusive Chaperone® Stretcher Exit System, gives control to the patient while helping reduce the risk of patient falls.

# Stryker Prime™ Series

stryker®



1105



1115



1125

## Medical

3800 E. Centre Ave.  
Portage, MI 49002 U.S.A.  
t: 269 329 2100 f: 269 329 2311  
toll free: 800 787 9537

www.stryker.com

### Standard Features

- Thermoformed ABS base with extra capacity storage tray
- 700 lb weight capacity
- 26-inch wide patient surface with 3-inch Enhanced Comfort Mattress
- Less than 21-inch low height
- Integrated utility tray
- Litter-mounted oxygen bottle holder
- Integrated transfer board
- Power-washable design (1105 and 1115 only)
- Four-wheel, center-locking brake system with dual end activators
- Dual side-mounted foot control hydraulics with uni-lower pedal
- Glideaway™ Siderails with integrated foot end push/pull handles
- 4 integrated IV receptacles
- +/- Trendelenburg
- 8-inch Omni Surface™ Casters with integrated wheel covers
- Pneumatic backrest/stationary foot
- Constant descent hydraulics
- Dual end brake/steer controls
- 2-year warranty

#### 1105: All Standard Features Plus

- Fifth Wheel steering

#### 1115: All Standard Features Plus

- Big Wheel steering

#### 1125: All Standard Features Plus C-US UL

- Zoom Motorized Drive System

### Optional Features

- 30-inch (76 cm) wide patient surface
- 4-inch (10 cm) Enhanced Comfort Mattress
- 4-inch (10 cm) or 5-inch (13 cm) Ultra Comfort Mattress
- Pioneer® Pressure Redistribution Mattress
- Four-sided brake and steer controls (1105 and 1115 only)
- Pop-up push handles (1105 and 1115 only)
- Knee Assist™ No-crank Hydraulic Knee Gatch
- Exclusive Lift Assist Backrest design with drop seat
- Integrated scale system
- Three-sided hydraulic controls (1105 and 1115 only)

### Specifications

	26-inch Litter	30-inch Litter
• <b>Model Number</b>	1105 (Prime with Fifth Wheel) 1115 (Prime with Big Wheel) 1125 (Prime with Zoom Drive)	1105 (Prime with Fifth Wheel) 1115 (Prime with Big Wheel) 1125 (Prime with Zoom Drive)
• <b>Overall Length</b>		
(1105 and 1115)	85" (216 cm)	85" (216 cm)
(1125)	86" (218 cm)	86" (218 cm)
• <b>Overall Width</b>		
Siderails Up	34" (86.4 cm)	38" (96.5 cm)
Siderails Down (1105 and 1115)	30.25" (76.8 cm)	30.5" (77.5 cm)
Siderails Down (1125)	30.25" (76.8 cm)	31" (78.7 cm)
• <b>Weight Capacity</b>	700 lbs (318 kg)	700 lbs (318 kg)
• <b>Height Range</b>		
High	34" (86.4 cm)	34" (86.4 cm)
Low (1105 and 1115)	20.75" (53 cm)	20.75" (53 cm)
Low (1125)	23" (58 cm)	23" (58 cm)
• <b>Litter Positioning</b>		
Backrest	0° – 90°	0° – 90°
Knee Gatch	0° – 40°	0° – 40°
Trend./Reverse Trend. (1105 and 1115)	±17°	±17°
Trend./Reverse Trend. (1125)	±16°	±16°
Maximum Angle	87.5°	87.5°
• <b>Patient Surface</b>	26" x 75.25" (66 x 191 cm)	30" x 75.25" (76 x 191 cm)
• <b>Siderails</b>	14" x 57.75" (36 x 147 cm)	14" x 57.75" (36 x 147 cm)
• <b>Caster Diameter</b>	8" (20 cm)	8" (20 cm)

### Warranty

- Two years parts, labor and travel

Stryker reserves the right to change specifications without notice.

### Electric Functionality

- Electric option (only 30 inches wide)
- Electric patient controls
- Foot end nursing controls
- Powered lift (1105 and 1115 only)
- Chaperone Stretcher Exit System
- 120V or 240V available

Copyright © 2010 Stryker  
PC/SCG 4M 0310, 09-SMDX-2842  
Mkt Lit-369-111709 Rev B