



76710-71M – Consists of :

- 76710 Transformer
- 23810 MacroView Otoloscope
- 11710 Ophthalmoscope

7670-12 – Mobile Stand

Model 23810 3.5v MacroView™ Otoloscope with Specula (Head Only)

- Twice the field of view.
- 30% greater magnification than a traditional otoscope.
- Ability to adjust focus to fine tune view.
- Tip grip for secure fastening of ear specula along with ejection feature for disposal.
- Fiber optics produce cool light with no reflections, no obstruction.



Model 11710 3.5v Standard Ophthalmoscope (Head Only)

- Halogen light for true tissue color and consistent, long-lasting illumination.
- 6 apertures for generalist use: micro, small, and large spot sizes, fixation target, slit aperture, and red-free filter.
- 28 focusing lenses with a range of -25 to +40 diopters.
- Rubber brow rest prevents scratching of eyeglasses.



Specifications

- Model No. **76720** 3.5v 767 Wall Transformer
(two-handle model with clock, instrument heads not included),
76710 3.5v 767 Wall Transformer
(two-handle model without clock, instrument heads not included),
76730 3.5v 767 Wall Transformer
Third Handle Module (instrument heads not included).

Power Cord:

Input: 108-125v., 60Hz. (0.5A @ 120V)


Output: Model Nos. 76720, 76710 — 3.00-3.90v, .700-1.05A at 120 VAC


Supplied with two coiled cords, two handles, and power cord.


Leakage current is less than 10 microamps from any exposed metal part.


Weight: Weight with mounting plate is 4.5 pounds (2.04kg). (NOTE: This is weight without instrument heads.)

Dimensions: H=4" (10.16cm) D=4" (10.16cm) W=12" (30.48cm)

 Listed to UL544

 Certified to CSA 22.2 No.125

 Certified

 The CE mark on this product indicates it has been tested to and conforms with the provisions noted within the 93/42/EEC Medical Device Directive.

Authorized European Representative Address:

European Regulatory Manager
Welch Alyn, Ltd.,
Kells Road, Navan,
County Meath,
Republic of Ireland
Tel. 353 46 28122
Fax 353 46 28536