



940 Cast Removal System

Cast Cutter Handpiece

PART #	DESCRIPTION
940	Cast Cutter

940 Cast Cutter Blades

PART #	DESCRIPTION
940-23	Ion Nitrided Blade (replaces standard 2 1/2" blade for use with synthetic or plaster casting material) (non-sterile) 1/pkg.
940-25	Stainless Steel Blade (replaces standard 2 1/2" blade) (non-sterile) 4/pkg.
940-35	Titanium Nitrided Blade (replaces standard 2 1/2" blade for use with synthetic casting material) (non-sterile) 2/pkg.
940-50	Cast Cutter Blade Tool (for blade removal on 940 Cast Cutter) 1/pkg.

Vacuums

PART #	DESCRIPTION
986	CastVac™ with 8 ft. hose and Mobile Stand (for use with 940, 840, and 848 Cast Cutters)
987	CastVac™ with 8 ft. hose and Wall Mount (for use with 940, 840, and 848 Cast Cutters)

Vacuum Accessories

PART #	DESCRIPTION
986-500	Replacement CastVac™ Bags
986-848	Adapter Cord and Hose (848 Cast Cutter to 986 CastVac™, 10 ft. cord, 8 ft. hose)
986-840	Adapter Hose (848 Cast Cutter to 986 CastVac™, 8 ft. hose)

Plaster Dispensers & Cast Room Organizers

PART #	DESCRIPTION
301	Hospital Plaster Dispenser, Stainless Steel (item 301 also includes the following accessories which may be purchased separately)
300-10-1	Disposable Buckets (non-sterile) 1/pkg.
300-30	Plastic Waste Basket 1/pkg.
300-40	Bucket Liners (non-sterile) 100/pkg.

General Accessories

PART #	DESCRIPTION
82	Small Cast Spreader
83	Large Cast Spreader
840-27	Cast Cutter Wall Bracket

Cast Walkers & Walking Heels

PART #	DESCRIPTION
93-100	Standard Adult Walking Heel 24/Pkg.
94-100	Child Size Walking Heel 24/Pkg.
98-100	Medium Size Walking Heel 24/Pkg.
99-100	SofTrac™ Sole Walking Heel 24/Pkg.
96	Toe Caps 12/Pkg.
97	Shower Shields 12/Pkg.

stryker
INSTRUMENTS



The Stryker 940 Cast Removal System

T H E N E W W O R L D S T A N D A R D

stryker
INSTRUMENTS

For more information, contact your local Stryker Sales Representative.

4100 East Milham Ave
Kalamazoo, MI 49001
(800)253-3210
(616)323-7700

940 Cast Cutter

Features/Benefits

- **Two Speed Control:**
 - **Standard Speed.** Very quiet, yet powerful enough to remove most casts.
 - **High Speed.** Increased speed for quickly removing even the thickest body jackets and hip spicas.
- **Versatile:** May be used with Vacuum or plugged directly into wall.
- **Quick Connect Vacuum Mount**
- **Impact Resistant Housing**
- **Traditional Ergonomic Design:** Ensures comfort and control while cutting.
- **Integral Vacuum Housing:** Reduces overall size while improving balance and mobility.



940 Cast Cutter shown actual size.

- **Efficient DC Motor:** With speed feedback, unparalleled power, less heat build up.

940 CastVac™

Features/Benefits

- **Detachable Cannister Lid:** Quick and easy emptying of dust collection filter.
- **Tool Bracket:** conveniently stores small cast spreader, scissors, screwdriver, and blade removal tool on back of stand.
- **Five Wheel Mobile Stand:** Provides stability and easy transportation.
- **Easy to operate Three Way Power Switch**
- **Twelve Foot Cord**



- **Quick Change Blade Mount**



- **Swivel Hose Mount:** Allows for easy movement of the 8ft. Vac Hose.
- **Quiet Vacuum:** Reduces noise level of the entire system.
- **Handpiece Control:** Vacuum automatically turns on when Cast Cutter is activated.
- **986 CastVac™ Specs:** 20" x 20" x 34.5" Tall

Cutting Accessories

940 Cast blades


Complete Cast Room Products

Fifty Years of Cast Room

LEADERSHIP



stryker®
INSTRUMENTS



Stryker Cast Room Products: "A continuing commitment to quality."



Complete Cast Room Products Ordering Information

CAST CUTTERS & ACCESSORIES

- 940 Cast Cutter
- 840 Standard 115V-Long Cord
- 848 Standard 115V-Short Cord (For use with all vacuums)
- 840-28 Carrying Bar
- 840-27 Wall Bracket
- 840-26 Wrench
- 840-4-2 Replacement Parts Catalog

VACUUMS & ACCESSORIES

- 864 OrthoVac®
- 986 CastVac™ with Mobile Stand
- 987 CastVac™ with 8ft hose and Wall Mount
- 988 Wall Mount only
- 989 Mobile Stand only
- 986-500 Filter Cartridge
- 940-886 Adapter Hose for use with 886 PlasterVac
- 986-840 Adapter Hose for use with 840 Cast Cutter
- 986-848 Adapter Hose for use with 848 Cast Cutter

BLADES

- 840-40-200 2 inch Stainless Steel Blade
- 840-40-250 2 1/2 inch Stainless Steel Blade
- 840-40-300 2 inch Titanium Nitrided Coated Blade
- 840-40-350 2 1/2 inch Titanium Nitrided Coated Blade
- 840-23 2 1/2 inch Ion Nitrided Blade
- 940-23 2 1/2 inch Ion Nitrided Blade
- 940-25 2 1/2 inch Stainless Steel Blade
- 940-35 2 1/2 inch Titanium Nitrided Blade



PLASTER DISPENSERS

- 301 Stainless Steel Dispenser (includes replacement parts)
- Replacement Parts**
- 300-10 Disposable Buckets (pkg. of 12)
- 300-20 Waxed Paper (300 ft roll)
- 300-30 Plastic Waste-Basket
- 300-40 Bucket Liners (pkg. of 12)

CAST SPREADERS

- 82 Cast Spreaders - Small
- 83 Cast Spreaders - Large

CAST WALKER & WALKING HEELS

- 93-100 Standard Adult Walking Heel (pkg. of 24)
- 98-100 Medium Size Walking Heel (pkg. of 24)
- 94-100 Child's Walking Heel (pkg. of 24)
- 99-100 SoffTrac® Sole Walking Heel (pkg. of 24)
- 97 Shower Shield (pkg. of 12)
- 96 Toe Caps (pkg. of 12)
- 806 Universal Cervical Collar

WARRANTY

Stryker Corporation products are warranted against defects in both material and workmanship for a period not to exceed one year from date of purchase.

For more information, contact your local Stryker sales representative.

Toll Free Order Department
1-800-253-3210
Fax: 1-800-999-3811

stryker®
INSTRUMENTS
4100 East Milham
Kalamazoo, MI 49001
1204