



SG-86-SS

MID-SIZE TABLES WITH SHELF

DIMENSIONS

SG-82-SS	20"W x 16"D x 34"H
SG-84-SS	30"W x 16"D x 34"H
SG-86-SS	33"W x 18"D x 34"H
SG-88-SS	36"W x 20"D x 34"H



SG-87-SS

MID-SIZE TABLES WITHOUT SHELF

DIMENSIONS

SG-83-SS	20"W x 16"D x 34"H
SG-85-SS	30"W x 16"D x 34"H
SG-87-SS	33"W x 18"D x 34"H
SG-89-SS	36"W x 20"D x 34"H



SG-92-SS

LARGE TABLES WITH SHELF

DIMENSIONS

SG-90-SS	48"W x 20"D x 34"H
SG-92-SS	48"W x 24"D x 34"H
SG-94-SS	60"W x 24"D x 34"H
SG-96-SS	72"W x 24"D x 34"H
SG-98-SS	36"W x 24"D x 34"H



SG-93-SS

LARGE TABLES WITHOUT SHELF

DIMENSIONS

SG-91-SS	48"W x 20"D x 34"H
SG-93-SS	48"W x 24"D x 34"H
SG-95-SS	60"W x 24"D x 34"H
SG-97-SS	72"W x 24"D x 34"H

FEATURES

- Stainless steel Steril-Gard construction eliminates full sub-tops where undesirable substances can collect
- Lightweight but strong
- Attractive design
- Easy to clean
- Soundproofed reinforcing channels for quiet use
- 3" rubber, double ball bearing casters for silent mobility
- Available in many sizes

OPTIONS

- 3 or 4 sided top and shelf rails available for all tables
- Stainless steel drawers, or tote boxes may be mounted under top.
- Table height can be modified

Glides/ Rubber Feet/ Casters

- #15 Glides
- #13 Rubber Feet
- #06 Locking Casters

Tote Boxes

Totes are constructed of impact-resistant, high-density polyethylene for durability

- CDS-TB-5" Tote Box (17¾"L x 14½"W x 5⅞"D)
- CDS-TB-10" Tote Box (17¾"L x 14½"W x 10"D)

Tote Box Hangers

- SG-TBH-1 16"D SG Tables
- SG-TBH-2 18"D SG Tables
- SG-TBH-3 20"D SG Tables
- SG-TBH-4 24"D SG Tables

IV Pole for SG Tables

- P-7201-1 available on lower shelf models only.

