

GO-CART DELUXE AND ACCESSORIES










Improve Point-of-Care Quality

- Secure storage with lockable stainless steel drawer
- Immediate access to all supplies during sample collection
- High quality stainless steel and aluminum construction
- Easily cleaned with standard disinfectants

See specifications tab for shipping restrictions.



DESCRIPTION	PREVIEW	UNIT	PRICE	EST. SHIP DATE	QTY
ML0379 Deluxe Package • Go-Cart Deluxe (as pictured) • 18.5"L x 22.5"W x 29"-45"H More Details		1/ea	\$2,060.00	Today	1
ML0380 Cart Only • Go-Cart Only • 18.5"L x 22.5"W x 29"-45"H More Details		1/ea	\$1,429.00	Today	1
ML0381 Accessory • Universal Bracket* (needed to mount containers) More Details		1/ea	\$152.00	3/5/2013	1
ML0382 Accessory • Locking Stainless Steel Drawer Assembly • OD: 17"W x 10"D x 5"H • ID: 15"W x 8"D x 3.5"H More Details		1/ea	\$327.00	Today	1
ML9166 Accessory • Stainless Steel Glove Box Holder More Details		1/ea	\$93.00	Today	1
ML0384 Accessory • 8qt Sharps Container with Mailbox Lid • 10.5"W x 7.25"D x 10"H More Details		20/pk	\$188.80	Today	1
ML0030 Accessory • 8qt Enviro Container • 10.5"L x 7.25"W x 10"H More Details		1/ea	\$13.20	Today	1
ML0031 Accessory • Mounting Bracket for 8qt Container • Non-locking • 10"W x 6.5"D x 7.5"H More Details		1/ea	\$20.90	Today	1

Details Specifications

A simple squeeze of a hand lever sets the Go-Cart height to a comfortable level when travelling, or to the height that positions supplies within easy reach in a patient room. The Go-Cart comfortably accommodates an 8qt sharps container, 8qt waste container, glove box, plus your phlebotomy tray with room to spare. Choose from a deluxe cart or build your own with separate accessories. To build your own Go-Cart, see the options under the "Related Products" tab.

*Universal Bracket is adjustable and compatible with a variety of wire style hangers for mounting containers.

- Go-Cart Deluxe
- Cart w/writing surface
 - ML1885 Phlebotomy Tray
 - Locking drawer
 - Enviro container
 - Sharps container