



Emergency

Standard Steel

6 Drawer Basic Emergency Cart: 30"

Beige, Breakaway Lock

Model: ABL-B-6



Standard features:

Standard 5" - 4 swivel, 2locking and 1 tracking heavy-duty casters and a built-in stabilizer bumper give Armstrong carts mobility and stability.

Overall Dimensions:

46 1/4" x 25" D x 33 1/2" W
(including bumpers,handles, casters)
117.5cm x 63.5cm x 85cm

Drawer Size:

17 1/6" D x 22 1/4"W
43.4cm Dx 56.5cm W

Drawer Configuration:

3"H (4), 6"H (1), 12"(1)
7.6cm H (4), 15.2cm H (1), 30.5cm (1)

Weight: 160 lbs / 72.5 kgs

Model ABL-B-6 (pictured)

Color: Beige

Note: Carts may require some assembly.

Frame color: Beige

LATEXfree

J.D. Honigberg International, Inc. 85 Revere Drive, Suite D, Northbrook, IL 60062 USA
Tel: 847-412-0200 • Fax: 847-412-0034 • medical@jdhintl.com



Cart Systems

- A-Smart Premier
- A-Smart Standard
 - Anesthesia / Procedure
 - ▶ Emergency
- Cart Accessories
- Utility/Stainless Steel
- Specialty

CPR and Training


Medical Equipment

Personal Protection

Anatomical Models



Printed Catalog
Order our full color printed catalog



Subscribe
Join our email list for news & updates

[Home](#) / [Cart Systems](#) / [A-Smart Standard](#) / [Emergency](#)




Standard Steel Two-Tone Emergency Carts

Standard 5" locking and swivel, heavy-duty casters and a built-in-stabilizer bumper give Armstrong carts mobility and stability, plus 18 colors to choose from (from top to bottom: SG - Sage Green, MG - Moss Grey, TC - Terra Cotta, SB - Slate Blue, L - Lavender, M - Mauve, Y - Yellow, CB - Taupe, T - Teal, B - Beige, AB - Crash Cart Blue, CR - Cranberry, BG - Blue Grey, HG - Hunter Green, **AA - Deep Blue**, AR - Red, S - Silver, CW - Cherry Woodgrain)



- Carts require some assembly.

Product Options ▾

	<p>Standard Steel Two-Tone 6-Drawer, Breakaway Lock, 30in. Product No: ABL-COLOR-6</p> <ul style="list-style-type: none"> ▫ Overall dimensions: 46 1/4"H x 25"D x 33 1/2"W (includes bumper/handles/casters) ▫ Drawer Size: 17 1/16"D x 22 1/4"W ▫ Drawer heights: 3"H (4), 6"H (1), 12"H (1) 	<p>Price: Call for Pricing</p>
	<p>Standard Steel Two-Tone 4-Drawer, Breakaway Lock, 24in. Product No: ABL-COLOR-4</p> <ul style="list-style-type: none"> ▫ Overall dimensions: 39 3/4"H x 25"D x 33 1/2"W (includes bumper/handles/casters) ▫ Drawer size: 17 1/16"D x 22 1/4"W ▫ Drawer heights: 3"H(2), 6"H(1), 12"H(1) 	<p>Price: Call for Pricing</p>
	<p>Standard Steel Two-Tone 5-Drawer, Breakaway Lock, 24in. Product No: ABL-COLOR-5</p> <ul style="list-style-type: none"> ▫ Overall dimensions: 39 3/4"H x 25"D x 33 1/2"W (includes bumper/handles/casters) ▫ Drawer size: Drawer size: 17 1/16"D x 22 1/4"W ▫ Drawer heights: 3"H (3), 6"H (1), 9"H (1) 	<p>Price: Call for Pricing</p>