

2.0 COMPONENT IDENTIFICATION

2.1 Spinal Surgery Top

The Spinal Surgery Top is 17" (43 cm) wide, 84" (213 cm) long, made from a carbon composite material and has a patient weight limit of 500 pounds (226 kg). The table top is labeled to indicate its patient weight limit. The Spinal Surgery Top includes the following items as standard components (see Figure 1):

5943 and 5943AP SPINAL SURGERY TOPS (Standard Components)		
Components	5943	5943AP
Small Hip Pad	2	2
Regular (Medium) Hip Pad	2	2
Large Hip Pad	4	2
Thigh Pad (Right/Left)	2	2
Chest Pad Assembly	1	1
Head Support Section (Face Plate)	1	1
Leg Support Section (Leg Boards)	2	2
Leg Sling	1	1
Buttock Strap	1	1
Safety Strap 60"	1	1
Safety Strap 90"	4	4
Articulating Arm Board	2	2

Figure 1: 5943/5943AP spinal surgery tops (standard components)




Figure 2: 5943 spinal surgery top with components
(The table top is shown coupled to the 5803 Advanced Control Modular Table System Base.)