



FS62L7BIMED

Maintain the lowest temperatures with energy efficiency in SUMMIT's innovative FS62L7BIMED, a slim all-freezer specially featured for use in medical, scientific, and laboratory environments.

The FS62L7BIMED has a slim width under 24 inches to fit in tight spaces. Designed for built-in installation, it can also function freestanding. The white cabinet includes a front-mounted lock for added security.

The manual defrost FS62L7BIMED operates several degrees colder than most freezers for safe long-term storage of food, vaccines, or pharmaceuticals. Temperatures can reach and maintain -25°C . Three removable pull-out basket drawers provide convenient storage solutions, while extra thick insulation limits energy waste and keeps the cold inside.

We include an audible temperature alarm complete with an external display mounted on the door. A green dot hospital grade plug is installed for added safety in high traffic areas.

With its slim size and reliable construction, the FS62L7BIMED works well in any setting in need of a trusted freeze. To create a complete undercounter set equipped for medical use, match the FS62L7BIMED

with SUMMIT's FF6L7BIMED all-refrigerator.

Highlights

- Alarm warns of temperature variations while providing a digital display of interior temperature
- Added safety with a hospital grade cord and plug
- Flexible design allows built-in or freestanding use
- Commercially approved for use in foodservice establishments
- Capable of operating at -25°C for long-term food preservation
- Perfectly sized to fit under counters and in tight 24 inch spaces
- Factory-installed lock provides security you can count on

Features

Commercially approved	ETL-S Approved to ANSI-NSF Standard 7 and CalCode Compliant
Slim undercounter dimensions	Less than 24 inches wide to fit tight spaces
Built-in capable	Make the best use of space by installing your appliance under the counter and flush with other cabinets
Fully finished white cabinet	Allows the freezer to be used freestanding
-25°C capable	Low temperature operation is ideal for vaccine, pharmaceuticals, and frozen foods
Factory installed lock	Keep contents secure with a key lock on the door
Alarm with temperature display	Temperature alarm warns of interior fluctuations while providing an external temperature readout on the door
Hospital grade cord and plug	Green dot plug for added safety
Manual defrost	No temperature fluctuations with an energy efficient manual defrost system
Three removable storage baskets	Conveniently organize your items while keeping the temperature stable when the door is opened
No internal fans	Static cooling system
Adjustable thermostat	Manage the interior temperature with ease and accuracy

One piece interior liner
100% CFC free

Enjoy easy clean-up with a seamless liner that won't hold a mess in hidden crevices
Environmentally friendly design without ozone-damaging chemicals

Specifications

Height:	33.5 up to 34.5"
Width:	23.63"
Depth:	23.5"
Capacity:	3.2 cu.ft.
Ext Depth with Handle:	25.13"
Interior Height 1:	24.25"
Interior Width 1:	17.75"
Interior Depth 1:	16.5"
Comp Step Height:	6.0"
Comp Step Width:	17.75"
Comp Step Depth:	6.0"
Height Adjustable:	Yes
Shipping Weight:	104.0 lbs.
Amps:	1.3
Canadian Electrical Safety:	ULC
US Electrical Safety:	UL
Reversible:	Factory Reversible
Door Swing:	RHD
Drawers:	Yes
Drawer Quantity:	3
Defrost Type:	Manual
Depth with door at 90°:	46.13"
Condensor Location:	Rear of Unit
Fan Type:	Exterior
Thermostat Type:	Dial
High Side PSI:	235.0
Level Legs Quantity:	4
Freon Type:	R134a
Low Side PSI:	70.0
Freon Oz:	2.7
Compressor Warranty:	5 Years
Parts/Labor Warranty:	1 Year
Sanitation:	ETL-S

Options

Cabinet Finishes	
Finishes	
White	
Door Finishes/Handles	
Finishes	
White	
Locks	
Standard Features	
L:	Single key lock
Upgrades	
KEYPAD:	Battery-powered keypad lock installed for keyless security
KEY:	Extra keys available for purchase
AC Keypad:	Keypad lock powered by appliance power cord for keyless security
Exterior Options	
Standard Features	
Alarm:	High/low temperature warning alarm with external temperature readout; meets most Health Dept. regulations
HG:	Hospital grade cord with green dot plug included
Upgrades	
Casters:	3" high NSF listed casters designed for institutional use, made from non-marring hard rubber in gray finish
TH:	Add a solar powered thermometer with external temperature readout
Logger:	USB Temperature Logger records interior temperature in 15 minute intervals to a USB flash drive for analysis on a PC or Mac
LHD:	Door reversal to left hand swing
Dolly:	Freestanding dolly with locking casters; ships fully assembled
SCD:	Protect interior contents and reduce energy waste with the Self Closing Door Option
Kickplate:	Add a stainless steel kickplate
DT:	Install a digital thermostat with external readout for precise temperature management
TT:	Install a battery-powered NIST traceable thermometer with external temperature readout; specify Fahrenheit or Celsius display upon ordering