



MEDASONICS®

Doppler Ultrasound Stethoscope

Features & Benefits

- **Advanced Audio Output**
- Offers superior volume control
- **Automatic Power-off**
- For extended battery life
- **Sensitec® Signal Processing**
- Provides pure, clear sound for greater accuracy in patient diagnosis
- **Blended Beam® Technology**
- Groundbreaking dual-crystal system designed for rapid signal acquisition
- **Accu-Rate™ Digital Display**
- Easy-to-read Accu-Rate Digital Heart Rate Display calculates beats per minute to save counting and calculating time
- **Excellent Value**
- A durable and elegant doppler at an affordable price
- **Flexible**
- Interchangeable fetal and vascular probes and optional waterproof fetal model provide the versatility you need in all clinical settings



CooperSurgical

Keeping You at the Forefront of Women's Health Care™

Doppler Ultrasound Stethoscope



Tria II Specifications

- Ultrasonic Frequency
- Fetal Doppler: 2.3 MHz
- Waterproof Fetal Doppler: 2.3 MHz
- Vascular Doppler: 5.0 MHz

Outputs

One 3.5 mm jack to connect to a headset or audio recording device

Controls

Volume control and push ON/OFF button

Power Source

One 9V lithium battery

Weight

8.75 oz. (245 gm)

Dimensions

5.25 x 2.6 x 1.75 in (13.3 x 6.6 x 4.5 cm)

Construction

Rugged ABS polycarbonate case

Service

CooperSurgical, Inc. maintains a full-service facility which can promptly repair products returned to the factory.

Warranty

Dopplers are covered by a 5-year warranty on parts, 1-year warranty on labor. Please call Customer Service for details.

Order Information

Description	Order Number	Frequency	Contents/Applications/Specifications
Tria II Fetal Doppler with Interchangeable Probe	T334	2.3 MHz	<u>Includes:</u> Doppler probe, integrated speaker unit, gel battery, soft carrying case and belt clip.
or with Dedicated Probe	T345	2.3 MHz	<u>Applications:</u> Detect fetal life early in pregnancy. Assess the rate and rhythm of the fetal heart.
Tria II Fetal Doppler with Accu-Rate Heart Rate Display and Interchangeable Probe	T334-D	2.3 MHz	<u>Includes:</u> Doppler probe, integrated speaker/unit, gel, battery, soft carrying case and belt clip.
or with Dedicated Probe	T345-D	2.3 MHz	<u>Applications:</u> Detect fetal life early in pregnancy. Assess the rate and rhythm of the fetal heart.
Tria II Vascular Doppler with Interchangeable Probe	T335	5.0 MHz	<u>Includes:</u> Doppler probe, integrated speaker/heart rate display unit, gel, battery, soft carrying case and belt clip.
or with Dedicated Probe	T346	5.0 MHz	<u>Applications:</u> Locate vessels and assess vessel patency and calculate ankle-brachial index.
Tria II Doppler with Interchangeable Fetal and Vascular Probes without heart rate display	T344	2.3 & 5.0 MHz	<u>Includes:</u> Doppler probes (2), integrated speaker unit (optional heart rate display), gel, battery, hard and soft carrying case and belt clip. <u>Applications:</u> See above Fetal and Vascular application descriptions.
with Accu-Rate Heart Display	T344-D		
Tria II Fetal Doppler with Waterproof Fetal Probe without heart rate display	T336	2.3 MHz	<u>Includes:</u> Waterproof Doppler probe, integrated speaker unit (optional heart rate display), gel, battery, soft carrying case and belt clip. <u>Applications:</u> Detect fetal life early in pregnancy. Assess the rate and rhythm of the fetal heart.
with Accu-Rate Heart Rate Display	T336-D		
Fetal Probe	P338	2.3 MHz	<u>Includes:</u> Fetal Doppler probe, gel.
Vascular Probe	P339	5.0 MHz	<u>Includes:</u> Vascular Doppler probe, gel.



Tria II Fetal Doppler



Tria II Fetal Doppler with Accu-Rate Heart Display



Fetal/Vascular Probe



Additional Accessories

