



**CHATTANOOGA MEDICAL SUPPLY, Inc.**

*The Professional's Choice*

Location: Hydrocollator Heating Units

**Model E-1 Stationary Heating Unit Special Pack**



Product Information

2102-2 - E-1 Stationary - Includes 3 Standard Size HotPacs, and 1 Cervical HotPac. 13"W x 8"D x 16"H (33 cm x 20 cm x 41 cm)  
(SMALLEST UNIT SHOWN)

- Power Supply: 110-120 Volts - 220-240 Volts
- Power Consumption: 1000 Watts
- Tank Capacity: 15.44 Liters
- Temperature Range: 160° - 165°F (71° - 74° C)
- Heating Up Time: 160°F (71°C) in 2 Hours
- Cool Down Time from 160°F (70°C): 1 Hour
- Safety Class 1: Type B
- Safety Tests: UL 544