

Instruction manual

MOTMed viva 2



12648 englisch

Technical Specifications, Symbols

Measurements and weight (basic model)

Height:	100 cm/40 inches
Width:	60 cm/23.6 inches
Depth:	56 cm/22 inches
Operating panel:	11.3 x 8.5 cm/4.5 x 3.4 inches
Weight: leg trainer	31 kg/68 lb
leg and arm trainer	38 kg/84 lb

Power requirements

Europe:	220 - 240 V ~/50 - 60 Hz, 0,56 A
USA, Canada:	110 - 120 V ~/50 - 60 Hz, 1,0 A
Japan:	100 - 110 V ~/50 - 60 Hz, 1,0 A
Other countries:	according to country specifications, consider the marking

Power consumption

Non-operating:	max. 7,5 VA
Operating:	max. 185 VA

Ambience conditions

Operation:	0° to +40°C/32 to 104°F
Storage:	-20° to +60°C/4 to 140°F
Humidity:	10 % - 80% relative, not condensing

System of protection: IPX0

Class of protection: I, Type B

**Class of medical products
according to MPG:** II a