

## LiteGait Jr.

### LiteGait® Jr. is ideal for larger size children.

The LiteGait Jr. model is perfect for a pediatric population, accommodating children up to 200 lbs. and 6 feet tall. The harness comfortably supports the child and the locking straps allow adjustments to be made for walking symmetry. The lighter design of the LiteGait Jr. allows a child to be safely supported while walking over ground or over a treadmill.



Gait therapy sessions using a treadmill at speeds as slow as 0.1 mph have proven highly effective with pediatric populations. In combination with the proper training program, LiteGait Jr. can help improve gait patterns, posture, balance, muscle strength, and overall endurance levels. The harness is adjustable and fits comfortably, providing for a very enjoyable walking therapy experience. LiteGait Jr. is being used successfully in homes, children centers and rehabilitation facilities around the world.

All LiteGait models are backed by a one year warranty\*, a responsive Service Department, and an extremely active Education Department conducting affordable on-site training and CEU accredited courses. When you purchase a LiteGait, you receive the service, the education, and the quality you and your clients deserve.

\* Policies may differ outside the US.

### Specifications<sup>1</sup>

	Jr 150	Jr 200
Max Patient Weight (lb)	150	200
Max Patient Height <sup>2</sup> (in)	66	72
Unit Max/Min Height (in)	72/41	78/41
Harness Fits Girths (in)	17-33	22-45
Base Inside Width (in)	27	30
Overall Length (in)	42	42

<sup>1</sup> Approximate and subject to change without notice.

<sup>2</sup> Max Patient Height is measured over 6 inch high treadmill.

[LiteGait photos](#)

[LiteGait Benefits](#)

[LiteGait Protocol](#)

### Product Pages

#### LiteGait

- [LGI200E](#)
- [LGI200P](#)
- [LGI260E/LGI360E](#)
- [LGI400E/LGI400ES](#)
- [LGI500E/LGI500ES](#)

#### LiteGait MX

- [LG100MX](#)
- [LG200MX](#)
- [LG300MX](#)

#### LiteGait Jr