



**PEDIGO**

*fine healthcare equipment*



Home

Products

**Stretchers**

About Pedigo

Find a Dealer

Find a Rep

Show Schedule

Contact Us

## Footstools and Step Stools

- ▶ Strong, durable, all welded construction
- ▶ Legs and handrails constructed of heavy 16-gauge 1" tube
- ▶ Slip resistant corrugated rubber top surface
- ▶ Rubber floor tips
- ▶ Tip resistant design

[P-10](#)

[P-1011-SS](#)

[P-1010-A-SS](#)

[P-11](#)

[P-1012-SS](#)

[P-18](#)

[P-1010-SS](#)

[P-10-A](#)

[P-18-A](#)



### Footstools

#### Chrome Models

##### P-10

- ▶ Dimensions: 12"D x 16"W x 8-1/4"H

##### P-11

- ▶ Dimensions: 12"D x 30"W x 8"H

#### Stainless Steel Models

##### P-1010-SS

- ▶ Dimensions: 12"D x 16"W x 8-1/4"H

##### P-1011-SS

- ▶ Dimensions: 12"D x 30"W x 5"H

##### P-1012-SS

- ▶ Dimensions: 12"D x 30"W x 8"H