



FOOTSTOOL

P-1010-SS Stainless Steel

- 12"D x 16"W x 8 1/4"H

P-1011-SS Stainless Steel

- 12"D x 30"W x 5"H

P-1012-SS Stainless Steel

- 12"D x 30"W x 8"H

P-10 Chrome

- 12"D x 16"W x 8 1/4"H

P-11 Chrome

- 12"D x 30"W x 8"H



FOOTSTOOL WITH HANDRAIL

P-1010-A-SS Stainless Steel

P-10-A Chrome

- 12"D x 16"W x 8 1/4"H
- 38 1/2"H with handrail

FEATURES

- Strong, all-welded steel
- Legs constructed of 1" O.D. 16-gauge tubular steel
- Handrail constructed of 1"O.D. tubular steel
- Tip resistant, rugged, and easy to maintain
- Corrugated, black, rubber top surface
- Rubber floor tips



DOUBLE STEP STOOL

P-18 Chrome

- 12"D x 16"W x 8 & 16"H

P-18-A Chrome

- same as P-18, with handrail
- 47"H with handrail



STACKABLE FOOTSTOOL

P-1015 Stainless Steel

- 14" D x 19" W x 6" H

FEATURES

- Extra large platform
- Slip resistant top
- Heavy gauge legs and handles
- Can stack safely for increased height and interlock to create larger platforms



TRANSPORT CART

P-1015-TC Stainless Steel

FEATURES

- Transport cart holds up to six P-1015 for storage or mobility
- 3" diameter swivel casters



ALL PEDIGO PRODUCTS ARE PROUDLY MADE IN THE USA