

6702/6702W

Oversized Wheelchair Scales®

Wheelchair Scales:

- [6002](#)
- [6102](#)
- [6702](#)
- [6702W](#)
- [6202](#)
- [6772](#)

- Wheelchair, chair and stand-on scales.
- Oversized platforms accommodate electric and extra-large wheelchairs.
- Automatic zero.
- Weight recall.
- Reweigh.
- Cordless.



Shown with optional height gauge and handrail. Wall-mounted height gauge also available.

SPECIFICATIONS



6702/6702W Oversized Wheelchair Scales

Platforms:	Low profile, mounted on large heavy-duty rubber wheels.
Model 6702:	28"w X 30"d X 2"h
Model 6702W:	28"w X 32"d X 2"h
Accuracy:	1/10 lb-100 gm.
Weighing Range:	1000 lb/454 kg.
Readout:	Digital LED display in pounds and kilograms, or kilograms only.
Power Source-Cordless:	6 D-size disposable batteries. Line cord power available.
Automatic Zero:	Returns to zero, ready for the next patient.
Weight Recall:	Recalls last weight.
Reweigh:	Recomputes weight while patient is on scale.
Audible Signal:	Indicates weight is in memory.
Low Battery Indicator:	Indicates battery needs replacing.
Weighing Mechanism:	Electronic load cells.
Successive Weighing:	Weigh a continuous stream of patients without re-zeroing.
Handrail:	(Accessory) Heavy-duty stainless steel.
Height Gauge:	(Accessory) Stainless steel, in/cm.
Printer:	(Accessory) Built-in, paper-tape, uses standard thermal paper.
Computer Capability:	RS-232 output (optional).

Our new oversized wheelchair scales incorporate all of the new features of our standard model including automatic zero, weight recall, reweigh and cordless operation. The digital readout electronically displays weight in pounds and kilos.

The oversized platforms are designed to weigh patients in electric and oversized wheelchairs. The platform height is still a mere 2" providing high stability and easy wheelchair access. With its access ramps folded up, the 6702 may also be used as a stand-on or chair scale.

Model 6702SP is available with a 35"W X 30"D X 2"H platform.

[Need more info?](#)

Click [here](#) to request a brochure, demo, prices or to have a Scale-Tronix representative contact you.

Or call:
1 800 873
2001