



# MAYO STANDS



## P-1065-SS MAYO STAND

- Removable, extra large, 20" x 25" stainless steel tray for easy cleaning
- Convenient foot-pedal operated lock adjusts height from 36" to 62"
- Base designed to move under low clearance equipment with maximum stability
- All welded stainless steel construction
- 2" double ball bearing swivel casters



## P-1066-SS MAYO STAND

- Same as P-1065-SS except for tray dimensions: 16-1/4" x 21-1/4" "

## P-1066-A-SS MAYO STAND NEW!

- Same as P-1066-SS except hand operated

## P-1068-SS MAYO STAND

- Same as P-1065-SS except for tray dimensions: 12-5/8" x 19-1/8" "

## P-1068-A-SS MAYO STAND NEW!

- Same as P-1068-SS except hand operated



## P-1069-SS MAYO STAND

- Same as P-1065 except with 16-1/4" x 21-1/4" tray and three leg base for increased mobility

## P-1069-A-SS MAYO STAND NEW!

- Same as P-1069-SS except hand operated



## P-66 MAYO STAND

- Removable 19-1/8"L x 12-5/8" W stainless steel tray on chrome plated frame
- Hand operated, friction lock adjusts height from 32" to 53"
- 2" double ball bearing swivel casters



## P-65 MEDICAL TRAY UNIT

- Removable 19-1/8"L x 12-5/8" W stainless steel tray on chrome plated frame
- Hand operated, friction lock adjusts height from 30" to 40"
- 2" double ball bearing swivel casters

4000 SE COLUMBIA WAY  
VANCOUVER WA 98661

800.822.3501

FAX 360.696.1700  
www.pedigo-usa.com