



**BASIN STANDS &
KICK BUCKETS**



P-78 chrome basin stand

single basin stand

P-78

Chrome
14"D x 14"W x 34"H

- One removable 13¼" diameter, 7 quart stainless steel basin
- Stainless steel construction

P-1078-SS

Stainless steel
14"D x 14"W x 34"H

FEATURES

- Elevated ring basin bracket
- Open frame design allows access from all sides
- 2" diameter double ball bearing swivel casters
- Mobile and easy to clean



P-1078-WIS-SS basin stand

single basin stand with shelf

P-1078-WIS-SS

Stainless steel
14"D x 14"W x 34"H

- One removable 13¼" diameter, 7 quart stainless steel basin
- Stainless steel construction

double basin stand with shelf

P-1079-WIS-SS

Stainless steel
14"D x 28"W x 34"H



P-79 chrome basin stand

double basin stand

P-79

Chrome
14"D x 28"W x 34"H

- Two removable 13¼" diameter, 7 quart stainless steel basin
- Stainless steel construction

double basin stand

P-1079-SS

Stainless steel
14"D x 28"W x 34"H



kick bucket

P-1020-SS

Stainless steel

Overall size
13½" diameter x 14"H

- Removable 12 quart, stainless steel basin

FEATURES

- Balanced, tip-resistant
- Electro-polished, stainless steel frame
- Perimeter bumper to protect walls
- 2" swivel casters for easy mobility



sponge receptacle

P-1025-SS

Stainless steel

Overall size
6¼" diameter x 8"H

- Removable 8½ quart, stainless steel basin

