

Blickman 7807SS-HB Baker Model, single-basin solution stand, Stainless-Steel

\$279.29

[Contact Us](#)

Model: M364-7807SS-HB
Shipping Weight: 18lbs
Manufactured by: Blickman

BLICKMAN MODEL 7807SS-HB BAKER SOLUTION STANDS have a single or double basin with a supporting frame connected by your choice of a bottom shelf or H-brace. The bracing and shelf, as well as all stainless welded construction, enhance the stand's strength and stability. Stainless steel basin has large 8.5 qt. capacity yet is unusually quiet when dropped into the circular ring. Rubber bumpers are placed strategically inside the ring to reduce noise from metal to metal contact. Travels quietly on two-inch swivel casters.

Dimensions	Feature	Basins
15"L x 15"W x 33"H	H-Brace	1

Model: 7807SS-HB Part: 0717807100



[larger image](#)

[larger image](#)