

**Stools and Chairs**

**S P E C I F I C A T I O N S**

**Classic Series Stools**

Seat Height:  
18" - 24.5" (274-278)  
19" - 25.5" (279)  
Seat Diameter: 16.25"  
Seat Cushion Thickness: 3.25"  
Back Cushion Dimensions:  
14" W x 7.5" H (275, 277)  
Caster Base Diameter: 21.5"

Weight Capacity:  
350 lbs. (274, 276)  
225 lbs. (275, 277, 278, 279)

**Value Series Stools**

Seat Height:  
18" - 24.5" (270-273)  
Seat Diameter: 16.25"  
Seat Cushion Thickness: 3.25"  
Back Cushion Dimensions:  
14" W x 7.5" H (271, 273)  
Caster Base Diameter: 23"  
Weight Capacity:  
300 lbs. (270, 272)  
225 lbs. (271, 273)

**Classic and Value Series Options**

Locking casters  
Soft rubber casters  
Glides  
6" Lab height extension  
Foot ring

**425/427 Air Lift Physician Stools**

Seat Height:  
17.5" - 24" (425)  
17" - 23" (427)  
Seat Cushion Dimensions:  
16.5" W x 17" D  
Seat Cushion Thickness: 3"  
Back Cushion Dimensions:  
13" W x 9.5" H  
(with back kit accessory)  
Caster Base Diameter: 23"  
Weight Capacity: 225 lbs.

**281 Blood Drawing Chair**

Seat Height: 20"  
Seat Cushion Dimensions:  
20" W x 16.5" D  
Seat Cushion Thickness: 2"  
Back Cushion Dimensions:  
20" W x 16.5" H  
Inside Arm Width: 21"  
Outside Arm Width: 32"  
Storage Drawer Capacity:  
8.5" W x 15.25" D x 4.5" H  
Weight Capacity: 400 lbs.

**680 Side Chair**

Seat Height: 17"  
Seat Cushion Dimensions:  
18" W x 19.5" D  
Seat Cushion Thickness: 2"  
Back Cushion Dimensions:  
18" W x 13.5" H  
Inside Arm Width: 17.5"  
(on model with arms)  
Outside Arm Width: 23.5"  
(on model with arms)  
Weight Capacity: 400 lbs.


**Upholstery Colors:**



Ritter: Navy, Dusty Blue, Moss, Clay, Pebble Grey, Shadow, Perfect Plum and Black



Midmark Premium Solid: Blueberry, Silver Sage, Tea Green, Fossil Grey, Cashmere Blue, Terra Cotta and Black



Midmark Premium Pattern: Belagio, Siena, Tuscany, Firenze, Capri and Treviso



Midmark Ultra-Premium: Ultra-Pacific Blue, Ultra-Spruce, Ultra-Olivine, Ultra-Arctic, Ultra-Sky, Ultra-Canyon and Ultra-Raven Wing

\*There will be an additional charge for Midmark Ultra-Premium.

**Midmark is an ISO 13485 and ISO 9001 Certified Company.**

For more information or a demonstration, contact your Ritter or Midmark Dealer, or contact Midmark Corporation, Versailles, Ohio. Call: 1-800-MIDMARK Fax: 1-800-365-8631 midmark.com

© Midmark Corporation 2005  
Products subject to improvement changes without notice.  
Litho in U.S.A. 007-0558-00 Rev. D1 (3/08)



*Classic and Value Series Backrest*



*Classic Series Glides*



*Value Series Glides*



*Classic and Value Series Locking Casters*



*Classic Series Soft Rubber Casters*



*Value Series Soft Rubber Casters*



*Classic Series 277 shown with optional 6" lab height extension and foot ring.*

**Ritter**<sup>®</sup>  
by MIDMARK