



Series 4024G Green-Line Treatment Table



- All Green-Line wood has no added urea-formaldehyde and exceeds CARB I and II standards.
- All urethane fabric is CAL 01350 compliant.
- PVC-free construction.
- H-Brace with laminate storage shelf..
- High pressure laminate legs, apron and stretchers.



Options

- 00* – Special Table Height
- 04* – Table Safety Strap
- 05* – Headrest, adj.
- 06* – Backrest, adj., 350 lbs. weight cap
- 07* – Nose Cutout
- 08* – Drawer (N/A with Opt 69)
- 09* – Armboard, adj
- 16* – Stirrups, adj
- 19* – Slide Shelf, laminate
- 23 – Paper Dispenser
- 25 – Paper Cutter
- 27 – Paper Dispenser & Cutter Combo
- 32 – Pillow, full size
- 35 – Pillow, small
- 69* – GasSpring Backrest
- 99L* – Left Hand Table Model


*Factory Installed

| MODEL | LENGTH | WIDTH | HEIGHT |
|-----------|--------|-------|--------|
| 4024G | 72" | 24" | 31" |
| 4024-030G | 72" | 30" | 31" |

COLORS:

Upholstery:

Green-Line PVC-free Urethane is Eco-friendly and CAL 01350 compliant



752 Indigo Blue 753 Flint Gray 754 Mushroom 755 Port 756 Beige

Laminate:

Durable - Scratch & Dent Resistant - Easy to disinfect



7909 Fusion Maple

Color accuracy may vary according to your monitor's settings. Color samples may be obtained via Customer Service.

COMPANY INFO:

About Us
Contact Us
Buy from Hausmann

SERVICES:

Literature
Media
Warranty Registration

HAUSMANN PRODUCTS:

www.hausmann.com
www.proteamttables.com
Overstock Garage

FOLLOW US:

