



Howard Instrument Table with Shelf

INSTRUMENT TABLES (HOWARD AND SAWYER MODELS). New and improved design offers complete 18-gauge, type-304, stainless-steel construction. Lightweight yet extremely durable. Sanitary closed top-hem prevents bacterial build-up. Sealed reinforcing channels provide structural rigidity. 1-1/4" diameter legs are fully welded for strength and stability. There is a 10" clearance between shelves on the Howard Models. The Sawyer Model has an H-brace in lieu of a second shelf. H-braces are fully welded around the legs with no crevices for dirt to collect. All tables travel on three-inch casters.

Options: Guard rail on three sides of top shelf.

Model Number	Part Number	Short Description	Dimensions	Ship Weight
783655	0117836000	Howard Instrument Table with Shelf	60"L 24"W 34"H	112



- Products
- What's New?
- Understanding Blickman Stainless Steel
- Custom Fabrication
- Care & Cleaning of Stainless Steel.
- Lifetime Limited Warranty

- Information
- Request
- FAQ
- Find a Dealer
- Technical Information
- Find a Sales Rep
- Trade Shows
- Become a Distributor
- Terms & Conditions
- Download a Catalog
- Employment
- Join our E-Mail List

Blickman Product Documentation

Click the documentation name below to open the document.

- [Entire Blickman Catalog](#)
PDF, 50193 KB
- [Blickman Tables](#)
PDF, 1055 KB, 5/16/2004