

[Home](#) > [Manufacturers](#) > [Chattanooga Group](#)

Adapta ADP-400 Hi-Lo Traction Table with Foot-switch & Casters - ADP40000212FC

Item #: ADP40000212FC
Brand: Adapta®
Manufacturer: Chattanooga Group
U/M: Each
Model: Four Section with Foot-switch & Casters
Color: Black

Price:
Available to order

Qty: 





Color sample:

Note: We are not responsible for differences in color nuances.

Product Description

The Adapta ADP-400 Hi-Lo Traction Table offers electric height adjustments from 21.5"-41.5" and has a 400 lb. lift capacity. The adjustable four-section table top has a head section that adjusts up or down +/-25°. The Pelvic/lumbar section separates smoothly, reducing friction during lumbar traction.

Functions:

- Four section
- Friction-free gliding lumbar section
- Turret mounted traction pedestal
- Head and Foot sections incline up to 90 grad
- 400 lb lifting capacity
- Electric High-Low
- Adjustable height range between 21.5" to 41.5"
- Foot control
- Optional hand control (item #66948 sold separately)
- Retractable casters
- Available in six colors
- Traction unit sold separately (item #4759)
- One year warranty

Technical Specifications:

- Mains Power: 120V, 50 Hz
- Current: 1.6 A
- Power Input: 200W
- Weight: 200 lbs
- Dimensions: 77"L x 27"W x 21.5" to 41.5"H
- Electrical Class: Class 1, Type B.

NOTE: Custom made product lead time - 2 weeks.

This product is listed in the following categories:
