

2000



MDS806300D

## Sets the Standard for Quality

Accessories			
<b>MDS85181U</b>	O <sub>2</sub> Holder	<b>MDS85196</b>	Anti-Fold Device
<b>MDS85183</b>	I.V. Pole	<b>MDS85197A</b>	Overhead Anti-Theft Device
<b>MDS85190</b>	O <sub>2</sub> / I.V. Combo	<b>MDS85189</b>	Anti-Tip Device

Chair	Arm Style	Leg Style	Seat W x D	HCPCS	Wt. Cap.
<b>MDS806100*</b>	Full Length Perm.	Permanent Footrest	18 x 16" (457 x 406mm)	K0001	300lb (136kg)
<b>MDS806150*</b>	Full Length Perm.	Detachable Footrest	18 x 16" (457 x 406mm)	K0001	300lb (136kg)
<b>MDS806200*</b>	Full Length Perm.	Detachable Elev. Legrest	18 x 16" (457 x 406mm)	K0001	300lb (136kg)
<b>MDS806250*</b>	Desk Length Remov.	Detachable Footrest	18 x 16" (457 x 406mm)	K0001	300lb (136kg)
<b>MDS806250DFLA**</b>	Full Length Perm.	Detachable Footrest	18 x 16" (457 x 406mm)	K0001	300lb (136kg)
<b>MDS806300DFLA**</b>	Full Length Perm.	Detachable Footrest	18 x 16" (457 x 406mm)	K0001	300lb (136kg)
<b>MDS806300*</b>	Desk Length Remov.	Detachable Elev. Legrest	18 x 16" (457 x 406mm)	K0001	300lb (136kg)
<b>MDS806400**</b>	Desk Length Remov.	Detachable Footrest	20 x 16" (508 x 406mm)	K0001 or K0002 E2201	300lb (136kg)
<b>MDS806450**</b>	Desk Length Remov.	Detachable Elev. Legrest	20 x 16" (508 x 406mm)	K0001 E2201 K0195	300lb (136kg)

\*To Order, Add The Letter Color Code Show Below Designating Color To The End Of The Item Number. \*\*Available Only In Standard Black Upholstery.

RBV (Ruby)
NVY (Navy)
D (Black)