



GRAHAM-FIELD HEALTH PRODUCTS

HEMOCARE PRODUCT CATALOG

Three-In-One Steel Commodes

- Versatile 3-in-1 design can be used as a bedside commode, raised toilet seat or toilet safety frame.
- Welded steel construction provides strength and stability.
- Durable, grey-coated steel frame is attractive and easy to clean.
- Back bar detaches to accommodate use as toilet safety frame.

- Sturdy plastic seat snaps on easily.
- Seat height readily adjusts in 1 inch increments.
- Comes complete with commode pail, cover and splash shield.
- Includes laminated color operating instructions.
- Limited lifetime warranty.



Platinum Collection

Increased seat depth provides added stability and comfort. Contoured waterfall armrests are designed to enhance user comfort and provide a good gripping surface. Designed to meet the requirements of HCPCS code: E0163.

7103 4/cs



Silver Collection

Contoured armrests enhance user comfort. Designed to meet the requirements of HCPCS code: E0163.

7103A-4 4/cs

Imperial Collection

Extra wide platform seat and base are comfortable for long periods of sitting. Commode opening is contoured and indented in the front to accommodate the needs of larger users. Removable, cushioned back bar enhances user comfort. Designed to meet the requirements of HCPCS code: E0168.

7446A-2 2/cs.



ITEM	Width Between Arms	Overall Width	Seat Width	Seat Depth	Seat Height Range	Weight Capacity	Replacement Tip
7103	18"	23"	14"	16 1/2"	16 3/4" - 27"	300 lb	LX9141
7103A-4	17 1/2"	24"	14"	15"	15 1/2" - 23 1/2"	300 lb	9290
7446A-2	23 1/2"	26"	23"	18"	16 1/2" - 23"	400 lb	9140

Three-In-One Aluminum Commodes

- 3-in-1 design is versatile and economical.
- Aluminum frame is lightweight and easy to transport.
- Contoured armrests enhance user comfort.
- Sturdy plastic seat snaps on easily.
- Comes complete with commode pail, cover and splash shield.
- Limited lifetime warranty
- Includes laminated color operating instructions.



Three-In-One Aluminum Commode with Back Bar and Casters

4" rear swivel locking casters. Designed to meet the requirements of HCPCS code: E0164.

2215A-2 2/cs

Three-In-One Aluminum Commode with Removable Back Bar

Designed to meet the requirements of HCPCS code: E0163.

2195A-4 4/cs

Three-In-One Aluminum Commode with Back Bar and 400 lb Maximum Weight Capacity

400 lb maximum weight capacity. Designed to meet the requirements of HCPCS code: E0168.

2190A 2/cs

Bariatric



ITEM	Width Between Arms	Overall Width	Seat Width	Seat Depth	Seat Height Range	Weight Capacity	Replacement Tip
2195A-4	18 1/2"	26"	13 1/2"	14 1/2"	17 1/2" - 21 1/2"	300 lb	LX9141
2215A-2	16 3/4"	21 1/4"	14"	16 3/4"	21 1/2"	300 lb	-
2190A	22 1/4"	27"	14"	16 3/4"	18" - 23"	400 lb	LX9141

Steel Drop Arm Three-In-One Commodes

- Arms drop independently below seat level for safe lateral transfers.
- Easy-to-release drop arm mechanism may be operated by patients lacking manual dexterity.
- Arms automatically lock in raised position.
- Removable back bar facilitates use as a bedside commode, raised toilet seat or toilet safety frame.

- Welded steel construction increases strength and stability.
- Cushioned armrests enhance user comfort.
- Comes complete with commode pail, cover and splash shield.
- Includes laminated color operating instructions.
- Limited lifetime warranty.



Platinum Collection

Watertight, heat-sealed padded seat is easy to clean and comfortable. Cushioned armrests enhance user comfort. Chrome-plated steel frame is scratch-resistant and easy to clean. Designed to meet the requirements of HCPCS code: E0165.

6437A 1 ea

Silver Collection

Durable, grey-coated steel frame is attractive and easy to clean. Cushioned armrests enhance user comfort. Sturdy plastic seat snaps on easily. Designed to meet the requirements of HCPCS code: E0165.

6433A-2 2/cs

Imperial Collection

Extra wide platform seat and base are comfortable for long periods of sitting. Commode opening is contoured and indented in the front to accommodate the needs of larger users. Cushioned back bar enhances user comfort. Designed to meet the requirements of HCPCS code: E0168. 400 or 600 lb maximum weight capacity.

6438A 400 lb 1 ea

6439A-2 600 lb 2/cs



ITEM	Width Between Arms	Overall Width	Seat Width	Seat Depth	Seat Height Range	Weight Capacity	Replacement Tip
6437A	18"	23"	14"	16 1/2"	19 1/2" - 23 1/2"	300 lb	9130
6433A-2	18"	23 3/4"	14"	16 3/4"	18 3/4" - 22 3/4"	300 lb	9290
6438A	25 1/2"	31"	23"	18 1/2"	19 1/2" - 23"	400 lb	9290
6439A-2	24 1/2"	32 3/4"	23"	18 1/2"	19 1/2" - 23 1/2"	600 lb	LX9141