



Standard Steel Red 6-Drawer, Breakaway Lock, 30in.

Product No: AR-6

Double-wall construction and ball bearing drawer slides add to the quality of all Armstrong carts.

- Overall: 46 1/4"H x 25"D x 33 1/2"W (Add side and top clearances)
- Drawer heights: 3" (4), 6", 12"
- Drawer size: 17 1/16"D x 22 1/4"W

Standard Steel Red 6-Drawer, Breakaway Lock, 30in.

Product No: AR-6

- Double Side-Wall Construction
- Push handles (2)
- TENTE Casters
- Locking casters (2)
- Swivel caster
- Tracking caster
- Stabilizing frame with bumper
- Ball Bearing Drawer Slides
- • Carts require some assembly.