



[Search Site](#)



[Shopping Cart](#)



[Mailing List](#)



[Request A Quote](#)

[Illuminators](#)

[Protective Apparel](#)

[Darkroom Items](#)

[Cassettes](#)

[Markers](#)

[Film](#)

[X-Ray Accessories](#)

[Wolf PACS](#)

[Digital Imaging](#)

[Vet Products](#)

[Contact](#)

[Home](#) » [Cassettes](#) >> [Cassette Carts, Racks, and Holders](#)

## Wolf Cassette Carts

### Heavy Duty Cassette Cart

Designed for rugged use, the Wolf Heavy Duty Cassette Cart is the perfect addition for the busy radiology department. Constructed of heavy gauge steel, with a tough washable white polyester powder finish, our heavy duty cassette cart rolls about easily on four sturdy rubber casters that lock in position. The heavy duty model has 3 levels divided into 5 compartments that hold any size cassette from 14" x 17" on down (maximum 125 cassettes).

Specifications: 17" D x 20" W. x 50" H. Shipping Weight: ~100 lbs

[Click here for Assembly Instructions](#)

### Wolf Cassette Carts

Constructed of gray anodized aluminum. The Wolf Cassette Carts carry up to 64 cassettes effortlessly. Built in dividers keep cassettes separated, and 1 inch tread strip prevents slippage. Choose from our bi-level or single level (not shown) models.

Bi-Level (64 Cassettes) 25514  
50" H x 32-1/2" L x 17" W. Shipping Weight: ~64 lbs.

Single-Level (32 Cassettes) 25513  
42" H x 32" L x 17" W. Shipping Weight: ~50 lbs.

### White Cassette Carts

All steel construction with powder coat finish, units have rubber mat inside each compartment and padded bar-style handles and scuff resistant tops.

Tall Cassette Carts (with casters) 80-138 -- Eight 6" wide compartments with 16 3/4" vertical clearance.  
24" L x 18" W x 45 1/2" H. Shipping Weight: ~120 lbs.

Tall Cassette Carts (without casters) 80-137 -- Stationary version of the Tall Cassette Cart.



Multiple Plate/Cassette Rack



Wolf Cassette Cart