

## Midland® Electric Parallel Bars



Base edges are tapered to accommodate wheelchairs up to 29" wide. Attractive stainless-steel handrails and powder-coated uprights provide a distinctive, durable finish that lasts for years. All Midland® Parallel Bars are backed by a Midland® exclusive, comprehensive 5-year warranty. 300 lb. weight capacity unless indicated otherwise. All motorized parallel bars are UL listed.

### Motorized Height and Width Parallel Bars

Handrail adjusts electrically for: height, 27" to 37" and width, 15" to 25-1/2".

Width between uprights 30-1/2".

Manual-override crank and removable abduction board included.

### Motorized Height Parallel Bars

Fits adults and children.

Handrail height adjusts electrically for adults, 27" to 37" and for children, 17" to 27". Brackets included to lower rails to child height.

Width adjusts manually from 16" to 24".

Manual-override crank and removable abduction board included.

### Bariatric Parallel Bars

Adjustable split handrails.

15'L with 45"W base.

Width adjusts 26" to 36-1/2 ". Height adjusts 27" to 37".

Cannot attach child rails.

500 lb. weight capacity.

### Pediatric Handrails (not shown)

Handrail height adjusts electrically, 17" to 27".

Width adjusts manually from 16" to 24".

1412E

 [Additional Images](#)

