



[Home](#) > [Patient Transfer](#) > [Transfer Boards](#) >

Static Treated - Transfer Board - White 3/16" x 22" x 72"



Our Price: \$110.00

Qty:

[add to cart](#)

Availability: Usually Ships within 1 Business Day

Product Code: 7799PTD-ST-WH

[★ Add to Wish List](#)



Description | [Instructions](#) | [Return Policy](#)

- The SafetySure® TransferEase™ Patient Mover makes transfers over uneven surfaces safe and simple.
- Semi-rigid polyethylene board is covered with an anti-static coating
- Anti-static surface protects patients from sticking to the surface while reducing static build-up.
- Ten handles for ease of positioning and can be used to hang the shifter for storage.
- Easily slides on a bed or stretcher pad.
- Radiolucent and can be left under the patient during procedures
- Prevents stress on the lower back of personnel when shifting patients
- All edges have been rounded
- 22" W x 72" H x 3/16" thick
- Item Weight: 9 lbs

