



Waste > Medical Waste

6144 Step-On Container



For hands-free sanitary waste management.

- Assists in complying with OSHA regulation; 29 CFR 1910.1030 "Occupational Exposure to Bloodborne Pathogens."
- Factory Mutual (FM) and California State Fire Marshal (CSFM) approved for fire safety.
- Tight-fitting, overlapping lid to help contain odors.
- Quiet, heavy-duty pedal for hands-free sanitary use.
- Tough, puncture-resistant, all-plastic construction.
- Optional rigid liners for 6143, 6144, 6145, 6146 containers provide a leak-proof base to help contain liquids.



CSFM



AVAILABLE COLORS

Order #	Color	Product UPC/ UCC Code
FG614400 WHT	WHT	086876057065 / 00086876057065
FG614400 RED	RED	086876057058 / 00086876057058
FG614400 BEIG	BEIG	086876057041 / 00086876057041



White
WHT



Beige
BEIG



Red
RED

SPECIFICATIONS

	U.S.	Metric
Length:	16 1/4 in	41.3 cm
Width:	15 3/4 in	40.0 cm
Height:	23 5/8 in	60.0 cm
Diameter:		
Volume Capacity [Nom]:	12 gal	45.4 L
Volume Capacity [Max]:		
Volume Capacity [Min]:		
Carton Length:		
Carton Width:		
Carton Height:		
Carton Cube:	3.73 ft3	0.11 m3
Ship Weight/Carton:	10 lb	4.54 kg
Pack Quantity:		1
Cartons Per Pallet:		18
Material: HDPE base/PP lid		
Process: Injection Molding		

ADDITIONAL INFORMATION:

Instructions: RCP_IS6143.pdf

Chemical Resistance Guide: chem.pdf

Products in Step-On Containers

Item #	Description	Length	Width	Height	Diameter	Volume Capacity
6142	Step-On Can with Rigid Liner			15 3/8 in	13 1/4 in	4 1/2 gal
6143	Step-On Container	16 1/4 in	15 3/4 in	17 1/8 in		8 gal
6144	Step-On Container	16 1/4 in	15 3/4 in	23 5/8 in		12 gal
6145	Step-On Container	19 3/4 in	16 1/8 in	26 1/2 in		18 gal