

HAND THERAPY TABLES



| | Length | Width | Height Range |
|---------------|--------|-------|--------------|
| 74-10S | 32" | 24" | 28"-43" |

Mobile Hand Therapy Table

- Spring loaded height adjustment with 2 hand control levers
- 1 1/8" Gray laminate top only
- Chrome base
- Removable casters for stationary application



| | Length | Depth | Height |
|---------------|--------|-------|---------|
| 74-110 | 45" | 20" | 31" |
| shelf | 20" | 20" | 30 1/2" |

Station without shelf unit

Hand Therapy Station with Shelf Unit

- 1 1/8" laminate top
- 3 shelf storage
- Natural finish only

| | Length | Width | Height |
|---------------|--------|-------|--------|
| 74-120 | 45" | 20" | 31" |

Hand Therapy Station

- Same as above without shelf unit

| | Length | Width | Height Range |
|----------------|--------|-------|--------------|
| 74-17G* | 36" | 36" | 24"-36" |

Soft Curve Gas Lift Table*

- 1 1/8" laminate top with sit-close cut-out on all sides
- Gas spring assist lever for easy height adjustment
- 30" base also available with optional casters (-081) (see page 74)

| | Length | Width | Height Range |
|----------------|--------|-------|--------------|
| 74-15K* | 36" | 19" | 23"-33" |

Hand Therapy Table*

- Height adjustable legs with self indexing spring loaded plunger
- Easy read, numbered height strip

| | Length | Width | Height Range |
|----------------|--------|-------|--------------|
| 74-14H* | 36" | 36" | 24"-36" |

Soft Curve Hydraulic Lift Table*

- 1 1/8" laminate top with sit-close cut-out on all sides
- Foot pump height adjustable base with swivel snug-top
- 28" chrome plated steel base



* Specify Maple, Gray or Natural top colors when ordering