

**PHILIPS IntelliVue MP60/70**

**APPROVED MOUNTING OPTIONS**

**FIXED HEIGHT MOUNTS**

*Provides lateral positioning, independent display/FMS swivel and display tilt*



**M-Series Support Arm Kit**  
AG-0018-21

Includes: 12" (30.5 cm) Pivot Arm, Mounting Plate for IntelliVue Table Top Adapter, Down Post for Flexible Module Server (FMS), and Channel Cover. Channel sold separately.



**M-Series Support Arm – Rail Kit**  
AG-0018-23

Includes: 12" (30.5 cm) Pivot Arm for 8-10 x 25-40 mm horizontal rail applications, Mounting Plate for IntelliVue Table Top Adapter, and Down Post for Flexible Module Server (FMS).

**VARIABLE HEIGHT MOUNTS**

*Provides variable height positioning - 17" (43.2 cm) Range, lateral positioning, independent display/FMS swivel, and display tilt*



**VHM™ Variable Height Support Arm Kit**  
AG-0018-25

Includes: VHM Arm, Down Post for Flexible Module Server (FMS), Mounting Plate for IntelliVue Table Top Adapter, and Channel Cover. Channel sold separately.



**VHM™ Variable Height Support Arm – Rail Kit**  
AG-0018-27

Includes: VHM Arm for 8-10 x 25-40 mm horizontal applications, Mounting Plate for IntelliVue Table Top Adapter, and Down Post for Flexible Module Server (FMS).

**KEYBOARD MOUNTING OPTION**



**MP60/70 Keyboard Mounting Bracket**  
AG-0018-86

Permits a keyboard to be added to MP 60/70 mounts. Can optionally accept a 9" (22.9 cm) Keyboard Plate WM-0023-39, or 17" (43.2 cm) Keyboard Plate WM-0010-60.

**ROLL STAND**

*Provides display tilt*



**Roll Stand Kit**  
AG-0018-60

Includes: Base/Casters/Post, Counterweight, Handle, Utility Basket, Mounting Plate for IntelliVue Table Top Adapter, and Flexible Module Server (FMS) mount.

**CAMLOCK MOUNT**

*Provides display swivel and tilt*



**Camlock Mount**  
AG-0018-83

Permits system to be mounted where Camlock rails are available. Includes: Camlock plate with Camlock lock lever/post, Mounting Plate for IntelliVue Table Top Adapter, and Flexible Module Server (FMS) mount.


**PEDESTAL MOUNT**

*Provides display swivel and tilt*



**"Duckfoot" Pedestal Mount**  
AG-0018-36

Freestanding system mount. Includes: Base/post, Mounting Plate for IntelliVue Table Top Adapter, and Flexible Module Server (FMS) mount.



DIMENSIONAL TOLERANCES AS FOLLOWS UNLESS OTHERWISE SPECIFIED: X.XXX: ±.010 X.XX: ±.025 X.X: ±.100		REV ECN# DATE APPRV. CMG DESIGNER: CHECKED BY:	VHC AND 12" (300 mm) SUPPORT ARM GCGX DWG NO. DATE: 8-13-08 REV. LEV: PT
GCGX PROPRIETARY THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION THAT IS GCGX CORPORATION INTELLECTUAL PROPERTY. DO NOT DISCLOSE TO OR DUPLICATE FOR OTHERS EXCEPT AS AUTHORIZED IN WRITING BY GCGX.			