

Model 505 - Temperature Management Unit



Features:

- Small size frees up more usable workspace
 - Light weight makes transport and set-up a breeze
 - Freestanding or easily attached to IV pole, bedrail or optional rolling cart
 - Built-in hour meter makes it easy to monitor usage for preventative maintenance
 - Unique snap-fit hose swivels at three points for easy blanket attachment and positioning
- Compatible with [241[®] blood/fluid warming set](#) - provides two warming therapies within one machine

Specifications

Dimensions:

13" h x 10" w x 11" d (33 x 25 x 28 cm)

Leakage Current:

Meets regulatory standards for leakage current

Weight:

13.6 lb (6.2 kg)

Filter:

High-efficiency 0.2 μ m filter

Operating Temperature:

High: 43° \pm 3°C (109.4° \pm 5.4°F)

Med: 38° \pm 3°C (100.4° \pm 5.4°F)

Low: 32° \pm 3°C (89.6° \pm 5.4°F)

Device Rating:

110-120 VAC, 60 Hz, 9.5 Amperes

220-240 VAC, 50 Hz, 4.5 Amperes

100 VAC, 50/60 Hz, 9.5 Amperes