

Medical Products and Services

Areas of Focus

Our Products

Services and Solutions

[Sign in](#) | [Register today](#)
 to order on-line

To order by phone:

888.444.5440


E-mail us if you have questions:

ambulatorycare@cardinal.com



Welch Allyn
PC-Based Resting ECG

- 2-way EHR connectivity to improve accuracy
- Get immediate access to patient data
- Reduce paperwork and errors
- Operates in thin-client environments
- Portable and featherweight


[Watch Video »](#)

[Surgery Center](#) / [Equipment](#) / [Carts](#) /

Uni-Cart® Emergency Carts

- Color coded for emergency code situations

(Waterloo Healthcare) These mobile carts store and transport emergency equipment needed for code situations. The electronically welded steel cabinets have prepierced sides for mounting optional accessories. The carts feature a wraparound rubber bumper, breakaway seal mechanism and nonconductive 5" - dia. swivel casters (one with brake).

The three-drawer carts have one 4"H, one 3"H and one 6 3/8"H drawer and an 11"H compartment. The five-drawer carts have one 4 "H, three 3"H and one 6 3/8"H drawer and an 11"H compartment.

(Three-Drawer Uni-Cart® Cart 61242-040 shown with Accessory Package EMER-PKG in photo. Accessory Package EMER-PKG available separately.)

[View more products by Waterloo Healthcare](#)

[Click to Enlarge](#)



Cat. No.	Mfr. No.	Color
3-Drawer Carts (26 1/2"W x 35"H x 18"D)		
61242-040	UCM-3PR	Red
UC-3PB	UCM-3PB	Blue
5-Drawer Carts (26 3/4"W x 41 3/8"H x 18 1/4"D)		
61242-055	UCM-5PR	Red
UC-5PB	UCM-5PB	Blue

[Three-drawer Uni-Cart® cart 61242-040 shown with accessory package EMER-PKG](#)

Cat. No.	Mfr. No.	Components
Accessory Package		
EMER-PKG	EMER-PKGM	2 utility hooks, IV pole, cardiac board and bracket, drawer trays, 100 plastic seals, rail and supports, collapsible side shelf, 10 drawer dividers, oxygen tank bracket