


SUPER BRIGHT SPOT™

Pre-focused spot available for OB/GYN's, urologists, family practitioners, hospital clinics and LDR's!

SUPER BRIGHT SPOT™ FEATURES:

- > 5185 footcandles (56,000 Lux)
- > 3700° K color temperature (natural coloring)
- > 20-watt halogen bulb with 1,000-hour average bulb life
- > 4"-5" spot size (much larger than fiber optic lights)
- > Dichroic reflector filters out infrared to keep the beam cool
- > High-quality, space-saving 4-caster base system – safe and flexible portable movement with locking feature in 2 of 4 casters
- > 115V model standard; 230V models available
- > 5-year warranty
- > UL 60601-1/60601-2-41 Certified
- > Made in the U.S.A.



| MODELS | | |
|---------|---------|----------------------------------------------------------------|
| 115V | 230V | |
| #SB40FL | #SB42FL | Super Bright Spot™ Floorstand (with casters) |
| #SB40W | #SB42W | Super Bright Spot™ with Extension Arm and Universal Wall Mount |

| ACCESSORIES | |
|-------------|--------------------------------|
| #0007000PK | Set of 4 ESX Replacement Bulbs |

NOTE: The optical data are nominal values based on measurements done according to the UL/IEC Performance standard at 120V at 18"

LIT-0002, Rev. 1



BURTON MEDICAL
 An ISO Certified Company
 21100 Lassen Street
 Chatsworth, CA 91311
 Phone: 800-444-9909, 818-701-8700
 Fax: 800-765-1770, 818-701-8725
 www.burtonmedical.com