

[Home](#) > [Diagnostic Tests](#) > [Pregnancy Tests](#) > [FSC-23-900-526](#)

Fisher HealthCare Sure-Vue Urine hCG Cassette Test Kit, 30/Pk

Item#: FSC-23-900-526



Features of the Fisher HealthCare Sure-Vue Urine hCG Test Kit:

- Greater than 99% accuracy.
- Sensitivity of 25mIU/mL.
- Results in three minutes.
- Visual, one-step test for the qualitative determination of human chorionic gonadotropin (hCG) in serum or urine to aid in early detection of pregnancy.
- Utilizes combination of monoclonal and polyclonal antibody reagent to selectively detect elevated levels of hCG.
- Positive specimens react with the specific antibody-hCG colored conjugate and form a colored line at the T (test) area.
- Absence of this colored line suggests a negative result.
- A colored line at the C (control) area will always appear regardless of the presence or absence of hCG.
- Addition of LH (300mIU/mL), FSH (1000mIU/mL) and TSH (1000 μ IU/mL) to negative or positive urine/serum showed no cross-reactivity with hCG.

- Store at room temperature.
- Includes: Test devices and Disposable specimen droppers.